


Diffraction Analysis

A demo of code for analysis of time series of diffraction patterns (DP).

The demo is only limited to fitting peaks' profiles to Gaussian forms for extracting peaks' intensity, sorting according to the position of the delay stage, normalizing and removing outliers (due to cosmic rays).

Files Scheme



Files Descriptions

read_coordinates.py – reads coordinates of peaks from a mouse click. Need to be done for each group of peaks before running the analysis (only once per experiment unless sample change/beam shift). The script generates files **coordinate*.txt**

experiment_info.txt – files with parameters of the analysis (numper of time points before excitation with the laser, rotation angle, sizes of integration windows, etc.)

vizualize.py – vizualization methonds

files_kit.py – methods for extracting data from external .txt files

IntensityWindow.py – class for integrating intensity within a window

scan_info.py – methods to operate with image series and to extract global (not peak-specific) information for each scan: time delay values, total intensity, camera background, etc

peak_info.py – methods to extract intensities of individual peaks, as well as sort and normalize them and remove outliers

analyze_script.py (MAIN SCRIPT) – extract intensities for each peak groups and record them into .csv file. Requires the files path, `experiment_info.txt`, `coordinate*.txt` for at least one peak.